

FEATURES THAT MAKE THE DIFFERENCE

High power 61-watt bi-polar radio frequency

High power bi-polar radio frequency achieves 50-600c temperature inside the dermis causes tissue coagulation and denaturation leads to neocollagenesis & neoelastogenesis.

Approved by US-FDA, K-FDA & CE

Dual Frequency 1 & 2Mhz

1 Mhz offers accurate frequency for Anti-aging treatments like face & neck wrinkles, fine lines, crow feet, sagging skin

2 Mhz offers precise frequency for acne scar, pores & stretch marks

Choice of Both Insulated & Non-Insulated Needle Sets

Insulated Micro-needle delivers wider RF emission only in dermis, while Non-insulated delivers RF effect from epidermis to dermis independent to skin types

Adjustable depth from 0.5mm to 3.5mm

Offers flexibility to treat different targeted areas

36 needles array with 0.3mm diameter

Optimum needle bundle can be used comfortably in small & large areas

One cartridge last the patient's entire treatment session

Offers unlimited shots leads to superior ROI

Variable Pulse duration of 100ms to 800ms

Provides personalize patient experience

Chevron-Circuit controlled RF Flow

Allows even delivery of the radio-frequency (RF Energy) prevents hot spots during treatment leads to decrease the pain and downtime

Robotic micro-motor system, advances the needle insertion more gently & comfortably leads painless treatment

Vivace is a first device to use Robotic micro-motor system, it is more comfortable than Hydraulic Motors or Solenoids

Red & Blue LEDs complement each treatment

Red stimulates collagen remodeling & Blue disrupts acne bacteria